

Product Datasheet

Recombinant Human C1QBP/HABP1 Protein, His Tag (orb1552810)

Description	Recombinant Human C1QBP/HABP1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	C1QBP
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human C1QBP(Leu74-Gln282) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q07021
Expiration Date	6 months from date of receipt.