

## Product Datasheet

### Recombinant Human BDH2 Protein, His Tag (orb1552851)

|                        |                                                                                                                            |
|------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human BDH2 Protein,<br>His Tag                                                                                 |
| <b>Species/Host</b>    | E. coli                                                                                                                    |
| <b>Reactivity</b>      | Human                                                                                                                      |
| <b>Endotoxins</b>      | Please contact our customer service<br>for more information.                                                               |
| <b>Target</b>          | BDH2                                                                                                                       |
| <b>Note</b>            | For research use only                                                                                                      |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE                                                                                             |
| <b>Source</b>          | This protein contains the human<br>BDH2(Met1-Leu245) was fused with<br>the N-terminal His Tag and expressed<br>in E. coli. |
| <b>Uniprot ID</b>      | <b>Q9BUT1</b>                                                                                                              |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                             |