

Product Datasheet

Recombinant Human BCL2L2/BCL-W Protein, His Tag (orb1552854)

Description	Recombinant Human BCL2L2/BCL-W Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	BCL2L2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human BCL2L2(Ala2-Thr172) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q92843
Expiration Date	6 months from date of receipt.