

Product Datasheet

Recombinant Human ATG10 Protein, T7 Tag & His Tag (orb1552892)

Description	Recombinant Human ATG10 Protein, T7 Tag & His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	ATG10
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human ATG10(Met1-Thr190) was fused with the N-terminal T7 Tag,C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9H0Y0
Expiration Date	6 months from date of receipt.