

## Product Datasheet

### Recombinant Human AKR1C3 Protein, His Tag (orb1552952)

<b>Description</b>	Recombinant Human AKR1C3 Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	AKR1C3
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human AKR1C3(Met1-Tyr323) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<b>P42330</b>
<b>Expiration Date</b>	6 months from date of receipt.