

Product Datasheet

Recombinant Mouse Acrp30/Adipoq Protein, His Tag (orb1552978)

Description	Recombinant Mouse Acrp30/Adipoq Protein, His Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Adipoq
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Adipoq(Glu18-Asn247) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q60994
Expiration Date	6 months from date of receipt.