

Product Datasheet

Recombinant Mouse 2b4/Slamf4/Cd244/Cd244a Protein, His Tag (orb1552993)

Description	Recombinant Mouse 2b4/Slamf4/Cd244/Cd244a Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Cd244a
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Cd244(Gln20-Asn221) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q07763
Expiration Date	6 months from date of receipt.