

Product Datasheet

Recombinant Human CKM/CKB Protein, Strep tag & His Tag (orb1553009)

Description	Recombinant Human CKM/CKB Protein, Strep tag & His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CKM/CKB
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CKM(Met1-Lys381) was fused with the C-terminal Strep Tag. This protein contains the human CKB(Met1-Lys381) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P06732, P12277
Expiration Date	6 months from date of receipt.