

Product Datasheet

Recombinant Human TPSAB1 Protein, His Tag (orb1553026)

Description	Recombinant Human TPSAB1 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	TPSAB1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TPSAB1(Met1-Pro275, Met1-Val133 & Ser135-Gly199 & Asp201-Pro275) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q15661
Expiration Date	6 months from date of receipt.