

Product Datasheet

Recombinant Human PADI4 Protein, TwinStrep Tag (orb1553034)

Description	Recombinant Human PADI4 Protein, TwinStrep Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PADI4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PADI4(Met1-Pro663) was fused with the N-terminal TwinStrep Tag and expressed in Baculovirus-Insect cells.
Uniprot ID	Q9UM07
Expiration Date	6 months from date of receipt.