

Product Datasheet

Recombinant Human PADI2 Protein, TwinStrep Tag (orb1553035)

Description	Recombinant Human PADI2 Protein, TwinStrep Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PADI2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PADI2(Met1-Pro665) was fused with the N-terminal TwinStrep Tag and expressed in Baculovirus-Insect cells.
Uniprot ID	Q9Y2J8
Expiration Date	6 months from date of receipt.