

Product Datasheet

Recombinant Human EIF4G/EIF4G1 Protein, His&Trx Tag (orb1553054)

Description	Recombinant Human EIF4G/EIF4G1 Protein, His&Trx Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	EIF4G1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human EIF4G1(Glu557-Lys646) was fused with the N-terminal His&Trx Tag and expressed in E. coli.
Uniprot ID	Q04637
Expiration Date	6 months from date of receipt.