

Product Datasheet

Recombinant Human CD264/TNFRSF10D Protein, His Tag (orb1553209)

Description	Recombinant Human CD264/TNFRSF10D Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	TNFRSF10D
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNFRSF10D(Met1-His211) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q9UBN6
Expiration Date	6 months from date of receipt.