

## Product Datasheet

# Recombinant Human CD179B/IGLL1 Protein, His Tag (orb1553260)

<b>Description</b>	Recombinant Human CD179B/IGLL1 Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	IGLL1
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human IGLL1(Met1-Ser213) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<b>P15814</b>
<b>Expiration Date</b>	6 months from date of receipt.