

Product Datasheet

Recombinant Human CD158K/KIR3DL2 Protein, His Tag (orb1553282)

Description	Recombinant Human CD158K/KIR3DL2 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	KIR3DL2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human KIR3DL2(Met1-His340) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P43630
Expiration Date	6 months from date of receipt.