

Product Datasheet

Recombinant Human CD127/IL7R Protein, His Tag (orb1553319)

Description	Recombinant Human CD127/IL7R Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	IL7R
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human IL7R(Met1-Asp239) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P16871
Expiration Date	6 months from date of receipt.