

Product Datasheet

Recombinant Human CD98/SLC3A2 Protein, His Tag (orb1553352)

Description	Recombinant Human CD98/SLC3A2 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	SLC3A2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SLC3A2(Arg206-Ala630) was fused with the N-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P08195
Expiration Date	6 months from date of receipt.