


Product Datasheet

AAV5 capsid protein VP1 antibody (orb155551)

Description	AAV5 capsid protein VP1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from Adeno-Associated Virus 5 capsid protein VP1 (21-120/724 aa)
Target	AAV5 capsid protein VP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80 kDa
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...