



Product Datasheet

ABCG2 antibody (orb155559)


Description	ABCG2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ABCG2. Please contact us for the exact immunogen sequence. The peptide is available as orb12938.
Target	ABCG2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	-100 µg (in 200 µl): 0.5 mg/ml-200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for one month and at -20°C for a longer time. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	72 kDa
Uniprot ID	Q80W57 , Q9UNQ0 , Q7TMS5
Entrez	9429
Dilution Range	IF/ICC: 1:100-500, IHC-P: 1:100-500, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



WB analysis of mouse brain (lane 1), rat...



IF image of rat brain tissue using anti-...



Immunohistochemical staining of paraffin...