



Product Datasheet

Adenovirus hexon protein antibody (orb155595)

Description	Adenovirus hexon protein antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from Human Adenovirus 4 hexon protein (651-750/936 aa)
Target	Adenovirus hexon protein
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105 kDa
Uniprot ID	Q9YVE5
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Brain Cell
U...



Flow
cytometric
analysis of
Brain Cell
U...