



Product Datasheet

RNF14 antibody (orb1556)


Description	RNF14 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human RNF14 (71-170/474 aa)
Target	RNF14
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	Q9UBS8
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat uterus...



Western blot analysis of mouse Placenta ...



Western blot analysis of mouse Cerebellu...