

Product Datasheet

AFAP1L2 antibody (orb155605)

Description	AFAP1L2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human AFAP1L2 (751-818/818 aa)
Target	AFAP1L2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90 kDa
Uniprot ID	Q8N4X5
Dilution Range	WB=1:500-2000, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.