

Product Datasheet

AKT (Phospho-Tyr473) antibody (orb155628)

Description	AKT (Phospho-Tyr473) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human AKT3 around the phosphorylation site of Tyr473 (FS(p-Y)SA)
Target	phospho-AKT3 (Tyr473)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57 kDa
Uniprot ID	P31749
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:25, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.