


Product Datasheet

Androgen Receptor (Phospho-Ser515) antibody (orb155659)

Description	Androgen Receptor (Phospho-Ser515) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Androgen Receptor around the phosphorylation site of Ser515 (YP(p-S)PT)
Target	phospho-Androgen Receptor (Ser515)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	99 kDa
CAS Number	b3
Uniprot ID	P10275
Dilution Range	IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of MCF7 cell using phospho-...