


Product Datasheet

AP2M1 (Phospho-Thr156) antibody (orb155675)

Description	AP2M1 (Phospho-Thr156) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human AP2M1 around the phosphorylation site of Thr156 (QV(p-T)GQ)
Target	Phospho-AP2M1 (Thr156)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q96CW1
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Sample:
Spleen
☐ mouse
Lysate at
40 ...