


Product Datasheet

APC1 (Phospho-Ser688) antibody (orb155682)

Description	APC1 (Phospho-Ser688) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human APC1 around the phosphorylation site of Ser688 (SL(p-S)PV)
Target	phospho-APC1 (Ser688)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	215 kDa
Uniprot ID	Q9H1A4
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat test...