

Product Datasheet

APC6 (Phospho-Thr580) antibody (orb155685)

Description	APC6 (Phospho-Thr580) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human APC6 around the phosphorylation site of Thr580 (II(p-S)PP)
Target	phospho-APC6 (Thr560)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72 kDa
Uniprot ID	Q13042
Expiration Date	12 months from date of receipt.