



Product Datasheet

APOBEC3B antibody (orb155694)


Description	APOBEC3B antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human APOBEC3B. Please contact us for the exact immunogen sequence. The peptide is available as orb12939.
Target	APOBEC3B
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	-100 µg (in 200 µl): 0.5 mg/ml-200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	46 kDa
Uniprot ID	F7EWS7, Q9UHL7
NCBI	001270411, 001257340
Entrez	9582
Dilution Range	IF/ICC: 1:100-600, IHC-P: 1:100-600, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



IF image of rat testis tissue using APOB...



IHC-P image of rat testis tissue using A...



Western blot analysis of rat lymph node ...