

Product Datasheet

Argonaute 3 antibody (orb155721)

Description	Argonaute 3 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human Argonaute 3/eIF2C3 (251-350/860 aa)
Target	Argonaute 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	97 kDa
Uniprot ID	Q9H9G7
Dilution Range	IF=1:100-500, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.