


Product Datasheet

ATP5G2 antibody (orb155761)

Description	ATP5G2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human ATP5G2 (67-141/141 aa)
Target	ATP5G2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7 kDa
Uniprot ID	Q06055
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human la...