

Product Datasheet

BCAR1 (Phospho-Tyr751) antibody (orb155803)

Description	BCAR1 (Phospho-Tyr751) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human BCAR1 around the phosphorylation site of Tyr751 (GQ(p-Y)EN)
Target	phospho-BCAR1 (Tyr751)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	93 kDa
Uniprot ID	P56945
Expiration Date	12 months from date of receipt.