

Product Datasheet

beta-Amyloid (Phospho-1-42) antibody (orb155836)

Description	beta-Amyloid (Phospho-1-42) antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human beta-Amyloid (672-713/770 or 1-42/42)
Target	beta-Amyloid (1-42)
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid: 0.01M PBS(pH7.4)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	affinity purified by Protein A
MW	4.4kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.