

Product Datasheet

BPK(Phospho-Tyr551) antibody (orb155859)

Description	BPK(Phospho-Tyr551) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Btk around the phosphorylation site of Tyr551 (DE(p-Y)TS)
Target	Phospho-BTK (Tyr551)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	76 kDa
Uniprot ID	Q06187
Expiration Date	12 months from date of receipt.