


Product Datasheet

Brefeldin A antibody (orb155873)

Description	Brefeldin A antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GBF1 (181-280/1859 aa)
Target	GBF1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	206 kDa
Uniprot ID	Q92538
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ki...