


Product Datasheet

BTN2A1 antibody (orb155883)

Description	BTN2A1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human BTN2A1 + BTN2A2 (21-120/527 aa)
Target	BTN2A1 + BTN2A2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57 kDa
Uniprot ID	Q7KYR7
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
mouse
spleno...