



Product Datasheet

C1orf141 antibody (orb155906)

Description	C1orf141 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human C1orf141 (101-200/400 aa)
Target	C1orf141
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	Q5JVB7
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human lung tissue.