


Product Datasheet

C1Q antibody (orb155963)

Description	C1Q antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human C1QA (61-160/245 aa)
Target	C1QA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	P02747
Dilution Range	WB=1:50-200, IHC-P=1:50-200, IHC-F=1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...