



Product Datasheet

C21orf33 antibody (orb156005)


Description	C21orf33 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human C21orf33 (101-200/268 aa)
Target	C21orf33
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	P30042
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg/Test, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat hear...



Flow cytometric analysis of hela Cell us...



Western blot analysis of human Hela lya...