


Product Datasheet

C2orf53 antibody (orb156034)

Description	C2orf53 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human C2orf53 (331-412/412 aa)
Target	C2orf53
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	Q53SZ7
Expiration Date	12 months from date of receipt.



Western blot analysis of Lovo Cell Lysat...