



Product Datasheet

C4orf52 antibody (orb156075)


Description	C4orf52 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human C4orf52 (101-168/168 aa)
Target	C4orf52
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7.7/18 kDa
Uniprot ID	Q8N5G0
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Heart lys...



Western blot analysis of Mouse Intestine...



Western blot analysis of mouse Muscle ly...