



Product Datasheet

C7orf64 antibody (orb156155)

Description	C7orf64 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human C7orf64 (1-100/367 aa)
Target	C7orf64
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	Q5RL73
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 Human Cell...



Western blot analysis of MCF-7 Human cel...