

Product Datasheet

C9orf84 antibody (orb156217)

Description	C9orf84 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human C9orf84 (1-100/1444 aa)
Target	C9orf84
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	165 kDa
Uniprot ID	Q5VXU9
Expiration Date	12 months from date of receipt.