

Product Datasheet

CACNA1B (Phospho-N type) antibody (orb156237)

Description	CACNA1B (Phospho-N type) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CACNA1B (N type) (101-200/2339 aa)
Target	CACNA1B (N type)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	257 kDa
Uniprot ID	Q00975
Expiration Date	12 months from date of receipt.