



Product Datasheet

CD368 Antibody (orb1562579)

Description	CD368 Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IP
Immunogen	CD368 ectodomain fused with human Fc
Target	CD368
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Uniprot ID	Q8WXI8
Entrez	338339
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocytes...



Flow cytometry surface staining pattern ...