

Product Datasheet

IL3RA Antibody (orb1562687)

Description	Purified Rabbit Anti Human IL3RA Polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 207-304aa of human IL3RA protein
Target	IL3RA
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Dilution Range	WB: 1:500-1:2,000. IHC: 1:100-1:200. ELISA: 1:20,000-1:80,000
Expiration Date	12 months from date of receipt.