



Product Datasheet

Anti-mCherry VHH antibody (orb1562797)

Description	Alpaca-derived anti-mCherry VHH single domain antibody (molecular weight: 14.3 kDa), expressed in E....
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % NaN ₃
Form/Appearance	Liquid
Concentration	1mg/mL
Storage	4°C; Do not freeze.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
MW	15 kDa
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Western blot analysis of 293T cell lysat...