



Product Datasheet

Anti-Her2 VHH antibody (orb1562800)

Description	Alpaca-derived anti-HER2 VHH single domain antibody (molecular weight: 13.7 kDa), expressed in E. coli...
Conjugation	Unconjugated
Tested Applications	IF, IHC
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % NaN ₃
Form/Appearance	Liquid
Concentration	1mg/mL
Storage	4°C; Do not freeze.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
MW	15 kDa
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunofluorescence analysis of paraforma...