



Product Datasheet

Anti-hTNFR1 VHH antibody (orb1562802)


Description	Alpaca-derived anti-TNFR1 VHH single domain antibody (molecular weight: 14.2 kDa), expressed in E. c...
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % NaN ₃
Form/Appearance	Liquid
Concentration	1mg/mL
Storage	4°C; Do not freeze.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
MW	15 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of human hepatocel...



Line graph illustrates about the Ag-Ab r...



Immunohistochemical analysis of paraffin...