

Product Datasheet

Anti-GFP VHH antibody, biotinylated (orb1562804)

Description	Alpaca-derived anti-GFP VHH single domain antibody, biotinylated (molecular weight: 15.82 kDa), expr...
Conjugation	Unconjugated
Tested Applications	IF
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % NaN ₃
Form/Appearance	Liquid
Concentration	1mg/mL
Storage	4°C; Do not freeze.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
MW	15 kDa
Expiration Date	12 months from date of receipt.