

Product Datasheet

Flag Tag (3A6) Mouse mAb (HRP) Antibody (orb1562806)

Description	Flag Tag (3A6) Mouse mAb (HRP) Antibody
Species/Host	Mouse
Reactivity	Other
Conjugation	HRP
Tested Applications	WB
Immunogen	DYKDDDDK synthetic peptide conjugated to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.

